


IL22 Antibody (R33818)

Catalog No.	Formulation	Size
R33818-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	53342
Localization	Secreted
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This IL22 antibody is available for research use only.



Western blot testing of human adipose lysate with IL22 antibody at 0.3ug/ml. Observed molecular weight: 19~25kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: Interleukin-17D, interleukin-27

Application Notes

Optimal dilution of the IL22 antibody should be determined by the researcher.

Immunogen

Amino acids EKDADSINSSIDK were used as the immunogen for this IL22 antibody.

Storage

Aliquot and store the IL22 antibody at -20oC.