


IL2RA Antibody / CD25 [clone RM418] (R20434)

Catalog No.	Formulation	Size
R20434-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM418
Purity	Protein A purified from animal origin-free supernatant
UniProt	P01589
Localization	Cell surface
Applications	Immunohistochemistry (FFPE) : 1:50-1:100
Limitations	This recombinant IL2RA antibody is available for research use only.



IHC staining of FFPE human tonsil tissue with recombinant IL2RA antibody at 1:200.



IHC staining of FFPE human Hodgkin lymphoma tissue with recombinant IL2RA antibody at 1:200.

Description

This antibody reacts to human CD25 (IL2Ra).

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant IL2RA antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the N-terminus of human Interleukin-2 receptor subunit alpha was used as the immunogen for the recombinant IL2RA antibody.

Storage

Store the recombinant IL2RA antibody at -20oC.