


IL7R Antibody (R33394)

Catalog No.	Formulation	Size
R33394-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	16197
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This IL7R antibody is available for research use only.



Western blot testing of NIH3T3 lysate with IL7R antibody at 2ug/ml. Predicted molecular weight: ~51/60-90kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Interleukin 7 receptor, CD127

Application Notes

Optimal dilution of the IL7R antibody should be determined by the researcher.

Immunogen

Amino acids DDHSFWCHSQLE were used as the immunogen for this IL7R antibody.

Storage

Aliquot and store the IL7R antibody at -20oC.