


IL7RA Antibody / Interleukin-7 receptor alpha / CD127 [clone 22I71] (RQ8889)

Catalog No.	Formulation	Size
RQ8889	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human, Mouse, Rat
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	22I71
Purity	Affinity chromatography
UniProt	P16871
Localization	Cell membrane, secreted
Applications	Western Blot : 1:500 Immunohistochemistry (FFPE) : 1:50
Limitations	This IL7RA antibody is available for research use only.



IHC staining of FFPE human tonsil tissue with IL7RA antibody, HRP-secondary and DAB substrate. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of 1) human HeLa, 2) human Jurkat, 3) human HepG2, 4) human MCF7, 5) rat spleen, 6) rat PC-12, 7) mouse spleen and 8) mouse RAW264.7 cell lysate with IL7RA antibody. Predicted molecular weight ~52kDa but can be observed at higher molecular weights due to glycosylation.