


ILK Antibody (R33551)


Catalog No.	Formulation	Size
R33551-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3611
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This ILK antibody is available for research use only.



Western blot testing of human 1) HEK293, 2) A431, 3) HeLa and 4) Jurkat cell lysate with ILK antibody at 1ug/ml. Expected molecular weight: 51-59 kDa.



Western blot testing of mouse NIH3T3 cell lysate with ILK antibody at 0.5ug/ml. Expected molecular weight: 51-59 kDa.

Description

Additional name(s) for this target protein: Integrin-linked kinase

Application Notes

Optimal dilution of the ILK antibody should be determined by the researcher.

Immunogen

Amino acids DMIVPILEKMQDK were used as the immunogen for this ILK antibody.

Storage

Aliquot and store the ILK antibody at -20oC.