


## Importin 8 Antibody (R35830)

Catalog No.	Formulation	Size
R35830-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10526
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Importin 8 antibody is available for research use only.



IHC testing of FFPE human tonsil with Importin 8 antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

## Description

Additional name(s) for this target protein: IPO8, RANBP8, RAN-binding protein 8

## Application Notes

Optimal dilution of the Importin 8 antibody should be determined by the researcher.

1. This Importin 8 antibody will detect isoforms 1 and 2.

## Immunogen

Amino acids SAFNFGTVPSNN were used as the immunogen for this Importin 8 antibody.

## Storage

Aliquot and store the Importin 8 antibody at -20oC.