

## INCA Antibody (R34818)

Catalog No.	Formulation	Size
R34818-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	440068
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This INCA antibody is available for research use only.



Western blot testing of human kidney lysate with INCA antibody at 0.5ug/ml. Predicted molecular weight: ~12kDa.

## Description

Additional name(s) for this target protein: Caspase recruitment domain family, member 17, CARD17, Caspase-1 inhibitor INCA

## Application Notes

Optimal dilution of the INCA antibody should be determined by the researcher.

## Immunogen

Amino acids NHLTTQDSQIVLPS were used as the immunogen for this INCA antibody.

## Storage

Aliquot and store the INCA antibody at -20oC.