

Integrin beta 5 Antibody (F43425)


Catalog No.	Formulation	Size
F43425-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43425-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P18084
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This Integrin beta 5 antibody is available for research use only.

250
130
95
72
55

Integrin beta 5 antibody western blot analysis in NCI-H460 lysate.



Integrin beta 5 antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human liver tissue.

Description

Integrin alpha-V/beta-5 is a receptor for fibronectin. It recognizes the sequence R-G-D in its ligand.

Application Notes

Titration of the Integrin beta 5 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 216-244 from the human protein was used as the immunogen for this Integrin beta 5 antibody.

Storage

Aliquot the Integrin beta 5 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.