


Interferon beta Antibody (F51190)


Catalog No.	Formulation	Size
F51190-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F51190-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Primate
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P01574
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50 Flow Cytometry : 1:10-1:50
Limitations	This Interferon beta antibody is available for research use only.



IHC analysis of FFPE human hepatocarcinoma with IFNB1 antibody



Flow cytometric analysis of HepG2 cells using IFNB1 antibody (blue) compared to a [negative control](#) (green). FITC-conjugated goat-anti-rabbit secondary Ab was used for the analysis.



Description

IFNB1 has antiviral, antibacterial and anticancer activities.

Application Notes

Titration of the Interferon beta antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 39-66 from the human protein was used as the immunogen for this Interferon beta antibody.

Storage

Aliquot the Interferon beta antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.