


Interferon gamma Antibody (R30566)

Catalog No.	Formulation	Size
R30566	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	P01579
Localization	Cytoplasmic, membranous, extracellular (secreted)
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Interferon gamma antibody is available for research use only.



Western blot testing of Interferon gamma antibody and Lane 1: recombinant human protein 10ng; 2: 5ng; 3: 2.5ng; 4: 1.25ng

Description

Interferon gamma is an inflammatory cytokine that has been implicated in the development of fibrosis in inflamed tissues. The production of IFN γ , which is under genetic control, can influence the development of fibrosis in lung allografts. It is also produced by natural killer(NK) cells and most prominently by CD8 cytotoxic T cells, and is vital for the control of microbial pathogens. Interferon gamma is believed to be crucial for host defence against many infections. Genetically determined variability in IFN γ and expression might be important for the development of tuberculosis. It activates human macrophage oxidative metabolism and antimicrobial activity. In addition to having antiviral activity, Interferon gamma has important immunoregulatory functions and plays an important role in the control of neointima proliferation.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the Interferon gamma antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

An amino acid sequence from the C-terminus of human IFNg (KRKRSQMLFRGRRASQ) was used as the immunogen for this Interferon gamma antibody.

Storage

After reconstitution, the Interferon gamma antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.