


IP10 Antibody / CXCL10 [clone AACI-3] (RQ5015)

Catalog No.	Formulation	Size
RQ5015	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AACI-3
Purity	Affinity purified
UniProt	P02778
Applications	Western Blot : 1:500-1:2000
Limitations	This IP10 antibody is available for research use only.



Western blot testing with IP10 antibody. Predicted molecular weight: 8-11 kDa.

Description

This antimicrobial gene encodes a chemokine of the CXC subfamily and ligand for the receptor CXCR3. Binding of this protein to CXCR3 results in pleiotropic effects, including stimulation of monocytes, natural killer and T-cell migration, and modulation of adhesion molecule expression. [RefSeq]

Application Notes

Optimal dilution of the IP10 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human IP10 / CXCL10 was used as the immunogen for the IP10 antibody.

Storage

Store the IP10 antibody at -20oC.