


IRF2 Antibody (R34678)


Catalog No.	Formulation	Size
R34678-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3660
Applications	Western Blot : 2-6ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This IRF2 antibody is available for research use only.



Western blot testing of Jurkat lysate with IRF2 antibody at 2ug/ml. Estimated molecular weight: ~39kDa but routinely observed at 39-50kDa.



Western blot of HEK293 lysate overexpressing IRF2 probed with IRF2 antibody (mock transfection in lane 1). Estimated/observed molecular weight: ~39/39-50kDa.

Description

Additional name(s) for this target protein: Interferon regulatory factor 2

Application Notes

Optimal dilution of the IRF2 antibody should be determined by the researcher.

Immunogen

Amino acids KTSDITQARVKSC were used as the immunogen for this IRF2 antibody.

Storage

Aliquot and store the IRF2 antibody at -20oC.