


IRF2BP1 Antibody (R35397)


Catalog No.	Formulation	Size
R35397-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26145
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This IRF2BP1 antibody is available for research use only.



Western blot testing of human liver lysate with IRF2BP1 antibody at 0.5ug/ml. Predicted molecular weight: ~62kDa but can be observed at up to 75kDa.



Western blot of nuclear NIH3T3 lysate (A) and of mouse (B) and rat (C) colon lysates at 1ug/ml. Predicted molecular weight: ~62kDa but can be observed at up to 75kDa.

Description

Additional name(s) for this target protein: Interferon regulatory factor 2 binding protein 1

Application Notes

Optimal dilution of the IRF2BP1 antibody should be determined by the researcher.

Immunogen

Amino acids QPPTQHRLVARNGE were used as the immunogen for this IRF2BP1 antibody.

Storage

Aliquot and store the IRF2BP1 antibody at -20oC.