


IRF4 Antibody (MUM1) (R31635)

Catalog No.	Formulation	Size
R31635	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA, 0.025% sodium azide
Gene ID	3662
Localization	Nuclear
Applications	Western Blot : 0.5-1ug/ml IHC (FFPE) : 0.5-1ug/ml
Limitations	This IRF4 antibody is available for research use only.



Western blot testing of IRF4 antibody and Lane 1: HeLa; 2: Jurkat lysate; Predicted molecular weight ~51 kDa, observed here at ~40 kDa.



Western blot testing of IRF4 antibody and recombinant human protein (0.5ng).



IHC-P: IRF4 antibody testing of human tonsil tissue.

Description

Interferon regulatory factor 4, also known as MUM1, is a protein that in humans is encoded by the IRF4 gene. It is located on 6p25.3. IRF4 is a transcription factor, and it is essential for the development of T helper-2 (Th2) cells, IL17 -producing Th17 cells, and IL9 -producing Th9 cells. In melanocytic cells, the gene may be regulated by MITF. IRF4 is a transcription factor that has been implicated in acute leukemia. This gene is strongly associated with pigmentation, sensitivity of skin to sun exposure, freckles, blue eyes, and brown hair color. Additionally, IRF4 inhibition is toxic to myeloma cell lines, regardless of transforming oncogenic mechanism.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the IRF4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Human partial recombinant protein (AA 272-451) was used as the immunogen for this IRF4 antibody.

Storage

After reconstitution, the IRF4 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.