


## IRF4 Antibody (R34209)


Catalog No.	Formulation	Size
R34209-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3662
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This IRF4 antibody is available for research use only.



Western blot testing of Jurkat lysate with IRF4 antibody at 0.5ug/ml. Predicted molecular weight: 51 kDa.



Western blot testing of human 1) Daudi and 2) Jurkat lysate with IRF4 antibody at 0.5ug/ml. Predicted molecular weight: 51 kDa.

## Description

Additional name(s) for this target protein: Interferon regulatory factor 4

## **Application Notes**

Optimal dilution of the IRF4 antibody should be determined by the researcher.

## **Immunogen**

Amino acids HISNPEDYHRSIR were used as the immunogen for this IRF4 antibody.

## **Storage**

Aliquot and store the IRF4 antibody at -20oC.