


IRF5 Antibody (F48241)


Catalog No.	Formulation	Size
F48241-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48241-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q13568
Applications	Western Blot : 1:1000 Immunofluorescence : 1:10-1:50 Flow Cytometry : 1:10-1:50
Limitations	This IRF5 antibody is available for research use only.



Western blot analysis of IRF5 antibody and Ramos lysate. Expected/observed molecular weight ~57 kDa.



Confocal immunofluorescent analysis of IRF5 antibody with A549 cells followed by Alexa Fluor 488-conjugated goat anti-rabbit IgG (green). DAPI was used as a nuclear counterstain (blue).



IRF5 antibody flow cytometry analysis of Ramos cells (green) compared to a negative control (blue). FITC-conjugated goat-anti-rabbit secondary Ab was used for the analysis.

Description

IRF5 is a member of the interferon regulatory factor (IRF) family, a group of transcription factors with diverse roles, including virus-mediated activation of interferon, and modulation of cell growth, differentiation, apoptosis, and immune system activity. Members of the IRF family are characterized by a conserved N-terminal DNA-binding domain containing tryptophan (W) repeats.

Application Notes

Titration of the IRF5 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 81-108 from the human protein was used as the immunogen for this IRF5 antibody.

Storage

Aliquot the IRF5 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.