


ITGA11 Antibody (R35168)

Catalog No.	Formulation	Size
R35168-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	319480
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ITGA11 antibody is available for research use only.



Western blot testing of mouse skeletal muscle lysate with ITGA11 antibody at 1ug/ml.
Predicted molecular weight: ~133kDa.

Description

Additional name(s) for this target protein: Integrin alpha 11

Application Notes

Optimal dilution of the ITGA11 antibody should be determined by the researcher.

1. Not suitable for ICC testing

Immunogen

Amino acids PTEEDLSHAPQRNH were used as the immunogen for this ITGA11 antibody.

Storage

Aliquot and store the ITGA11 antibody at -20oC.