


ITK Antibody / IL2-inducible T-cell kinase (R34415)

Catalog No.	Formulation	Size
R34415-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3702
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This ITK antibody is available for research use only.



Western blot testing of rat thymus lysate with ITK antibody at 0.3ug/ml. Predicted molecular weight: ~72kDa.

Description

Additional name(s) for this target protein: Tyrosine-protein kinase ITK/TSK

Application Notes

Optimal dilution of the ITK antibody should be determined by the researcher.

Immunogen

Amino acids KERPEDRPAFSR were used as the immunogen for this ITK antibody.

Storage

Aliquot and store the ITK antibody at -20oC.