


Jagged 1 Antibody (R36227)

Catalog No.	Formulation	Size
R36227-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	182
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Jagged 1 antibody is available for research use only.



IHC testing of FFPE human kidney with Jagged 1 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: JAG1

Application Notes

Optimal dilution of the Jagged 1 antibody should be determined by the researcher.

Immunogen

Amino acids TNKQDNRDLESAQS were used as the immunogen for this Jagged 1 antibody.

Storage

Aliquot and store the Jagged 1 antibody at -20oC.