


JAM-A Antibody / F11R (R35400)

Catalog No.	Formulation	Size
R35400-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	50848
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This JAM-A antibody is available for research use only.



Western blot testing of human peripheral blood mononucleocyte lysate with JAM-A antibody at 0.3ug/ml. Expected molecular weight: 35~43kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: Junctional adhesion molecule A, F11 receptor, CD321; F11R

Application Notes

Optimal dilution of the JAM-A antibody should be determined by the researcher.

Immunogen

Amino acids QPSARSEGEFKQTS were used as the immunogen for this JAM-A antibody. This sequence is common to

isoforms 1 and 2.

Storage

Aliquot and store the JAM-A antibody at -20oC.

References (1)