


## JNK2 alpha Antibody (R33880)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33880-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                               |
|----------------------|---------------------------------------------------------------|
| Availability         | 1-3 business days                                             |
| Species Reactivity   | Human                                                         |
| Predicted Reactivity | Dog, Mouse, Rat                                               |
| Format               | Antigen affinity purified                                     |
| Clonality            | Polyclonal (goat origin)                                      |
| Isotype              | Goat Ig                                                       |
| Purity               | Antigen affinity                                              |
| Gene ID              | 5601                                                          |
| Applications         | Western Blot : 0.03-0.1ug/ml<br>ELISA (peptide) LOD : 1:16000 |
| Limitations          | This JNK2 alpha antibody is available for research use only.  |



Western blot testing of HeLa lysate with JNK2 alpha antibody at 0.03ug/ml. Predicted molecular weight: ~48kDa.

## Description

Additional name(s) for this target protein: C-Jun kinase 2 antibody; MAPK9

## Application Notes

Optimal dilution of the JNK2 alpha antibody should be determined by the researcher.

1. This JUNK2 alpha antibody will detect isoforms alpha1, alpha2. and gamma.

## Immunogen

Amino acids ELVKGCVIFQGDH were used as the immunogen for this JNK2 alpha antibody.

## Storage

Aliquot and store the JNK2 alpha antibody at -20oC.