


JUN Antibody / c-Jun [clone CAO-10] (RQ5122)

Catalog No.	Formulation	Size
RQ5122	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	CAO-10
Purity	Affinity purified
UniProt	P05412
Applications	Western Blot : 1:500-1:2000
Limitations	This JUN antibody is available for research use only.



Western blot testing of mouse NIH3T3 cell lysate with JUN antibody. Predicted molecular weight ~36 kDa, commonly observed at 36-43 kDa.

Description

This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression.[RefSeq]

Application Notes

Optimal dilution of the JUN antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human c-Jun / JUN was used as the immunogen for the JUN antibody.

Storage

Store the JUN antibody at -20oC.