


## JUNB Antibody (R36209)

Catalog No.	Formulation	Size
R36209-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3726
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This JUNB antibody is available for research use only.



Western blot testing of human breast cancer lysate with JUNB antibody at 0.3ug/ml.  
Expected molecular weight: 36-39 kDa (non-phosphorylated), 40-45 kDa (phosphorylated).

## Description

Additional name(s) for this target protein: Jun B proto-oncogene, AP-1, activator protein 1

## Application Notes

Optimal dilution of the JUNB antibody should be determined by the researcher.

1. This JUNB antibody may cross react with JUN/c-JUN protein.

## Immunogen

Amino acids TLKAENAGLSSTAG were used as the immunogen for this JUNB antibody.

## Storage

Aliquot and store the JUNB antibody at -20oC.