

JUND Antibody (F40151)

Catalog No.	Formulation	Size
F40151-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40151-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P17535
Applications	Western Blot : 1:1000
Limitations	This JUND antibody is available for research use only.



Western blot analysis of JUND antibody and Jurkat lysate. Predicted molecular weight: ~39/34kDa (JUND-L/JUND-S).

Description

JUND is a member of the JUN family, and a functional component of the AP1 transcription factor complex. It has been proposed to protect cells from p53-dependent senescence and apoptosis.

Application Notes

Titration of the JUND antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 233-262 from the human protein was used as the immunogen for this JUND antibody.

Storage

Aliquot the JUND antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.

References (1)