


JUND Antibody (R33312)

Catalog No.	Formulation	Size
R33312-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3727
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This JUND antibody is available for research use only.



Western blot testing of NIH3T3 nuclear lysate with JUND antibody at 1ug/ml. Predicted molecular weight: ~39/34kDa (JUND-L/JUND-S). (2)

Description

Additional name(s) for this target protein: Transcription factor jun-D

Application Notes

Optimal dilution of the JUND antibody should be determined by the researcher.

1. This JUND antibody is not expected to cross-react with Jun.

2. The short form of the protein lacks the N-terminal 43 amino acids.

Immunogen

Amino acids CQLLPQHQPAY were used as the immunogen for this JUND antibody.

Storage

Aliquot and store the JUND antibody at -20oC.

References (1)