


Kallikrein 6 Antibody (R34403)

Catalog No.	Formulation	Size
R34403-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	19144
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Kallikrein 6 antibody is available for research use only.



Western blot testing of mouse liver lysate with Kallikrein 6 antibody at 2ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight: ~27kDa.

Description

Additional name(s) for this target protein: Kallikrein related-peptidase 6, Klk6, kallikrein 7, Klk7

Application Notes

Optimal dilution of the Kallikrein 6 antibody should be determined by the researcher.

Immunogen

Amino acids KEKPGVYTDVCTHIR were used as the immunogen for this Kallikrein 6 antibody.

Storage

Aliquot and store the Kallikrein 6 antibody at -20oC.