


## KCNRG Antibody (F40014)


| Catalog No.   | Formulation                                | Size    |
|---------------|--------------------------------------------|---------|
| F40014-0.4ML  | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.4 ml  |
| F40014-0.08ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.08 ml |

[Bulk quote request](#)

|                           |                                                         |
|---------------------------|---------------------------------------------------------|
| <b>Availability</b>       | 1-3 business days                                       |
| <b>Species Reactivity</b> | Mouse                                                   |
| <b>Format</b>             | Antigen affinity purified                               |
| <b>Clonality</b>          | Polyclonal (rabbit origin)                              |
| <b>Isotype</b>            | Rabbit Ig                                               |
| <b>Purity</b>             | Antigen affinity                                        |
| <b>UniProt</b>            | Q8N5I3                                                  |
| <b>Applications</b>       | IHC (Paraffin) : 1:10-1:50<br>Western Blot : 1:1000     |
| <b>Limitations</b>        | This KCNRG antibody is available for research use only. |



KCNRG antibody analysis in formalin fixed and paraffin embedded human brain tissue.



KCNRG antibody western blot analysis in mouse liver, lung tissue lysate

## Description

This gene encodes a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This gene is located on chromosome 13 within the minimal deletion region for B-cell chronic lymphocytic leukemia. Multiple alternatively spliced transcript variants have been found for this gene.

## Application Notes

Titration of the KCNRG antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 175-202 from the human protein was used as the immunogen for this KCNRG antibody.

## Storage

Aliquot the KCNRG antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.