


KIAA0319 Antibody (R35274)

Catalog No.	Formulation	Size
R35274-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	210108
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This KIAA0319 antibody is available for research use only.



Western blot testing of mouse brain lysate with KIAA0319 antibody at 0.5ug/ml. Both the expected ~150kDa band and additional ~25kDa band are blocked by immunizing peptide.

Description

Additional name(s) for this target protein: D130043K22Rik, 4930451E12Rik

Application Notes

Optimal dilution of the KIAA0319 antibody should be determined by the researcher.

Immunogen

Amino acids QGKIKQENKPTLH were used as the immunogen for this mouse KIAA0319 antibody.

Storage

Aliquot and store the KIAA0319 antibody at -20oC.