


KIF3A Antibody (R32364)

Catalog No.	Formulation	Size
R32364	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide
UniProt	Q9Y496
Applications	Western Blot : 0.1-0.5ug/ml
Limitations	This KIF3A antibody is available for research use only.



Western blot testing of 1) rat brain, 2) rat testis, 3) mouse brain, 4) mouse testis, 5) human MCF7 lysate with KIF3A antibody. Expected/observed molecular weight ~80 kDa.

Description

Kinesin-like protein KIF3A is a protein that in humans is encoded by the KIF3A gene. KIF3A is one subunit of the heterotrimeric motor protein, kinesin-2, that was initially isolated from sea urchin egg/embryo cytosol using microtubule affinity purification. This motor consists of two kinesin-related subunits (called KIF3A and KIF3B or 3C in vertebrates) and an associated protein (KAP3), and it transports protein complexes, nucleic acids and organelles towards the plus ends of microtubule tracks within cells. Work done in a broad range of eukaryotic cells has revealed that heterotrimeric kinesin-2 is the primary motor protein driving the intra-flagellar transport of tubulins and other axonemal building blocks from the base of the ciliary/flagellar axoneme to their site of assembly at the distal tips. This process is required for cilium assembly/maintenance and cilium-based signalling which play key roles in various cell and developmental processes. For

example, in vertebrate embryos, kinesin-2 function is required for cilia-dependent nodal flow and the development of left-right asymmetry.

Application Notes

Optimal dilution of the KIF3A antibody should be determined by the researcher.

Immunogen

Amino acids 485-699 of human KIF3A were used as the immunogen for the KIF3A antibody.

Storage

After reconstitution, the KIF3A antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.