


Klotho Antibody (R34609)


Catalog No.	Formulation	Size
R34609-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9365
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This Klotho antibody is available for research use only.



IHC testing of FFPE human breast with Klotho antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.



IHC staining of FFPE human placenta with Klotho antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: KL

Application Notes

Optimal dilution of the Klotho antibody should be determined by the researcher.

Immunogen

Amino acids KRDDAKYMYYLKK were used as the immunogen for this Klotho antibody.

Storage

Aliquot and store the Klotho antibody at -20oC.