


KMO Antibody / Kynurenine 3 monooxygenase (R36453)

Catalog No.	Formulation	Size
R36453-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	O15229
Gene ID	8564
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This KMO antibody is available for research use only.



Western blot testing of human placenta lysate with KMO antibody at 0.03ug/ml. Predicted molecular weight: ~56 kDa. Predicted molecular weight: 52-56 kDa (multiple isoforms).

Description

Additional name(s) for this target protein: Kynurenine 3-hydroxylase

Application Notes

Optimal dilution of the KMO antibody should be determined by the researcher.

Immunogen

Amino acids NTTCTPAKAVDSLE were used as the immunogen for this KMO antibody.

Storage

Aliquot and store the KMO antibody at -20oC.