


Kynurenine 3 monooxygenase Antibody / KMO (R36454)

Catalog No.	Formulation	Size
R36454-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	O15229
Gene ID	8564
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Kynurenine 3 monooxygenase antibody is available for research use only.



Western blot testing of human placenta lysate with Kynurenine 3 monooxygenase antibody at 1ug/ml. Predicted molecular weight: 52-56 kDa (multiple isoforms).

Description

Additional name(s) for this target protein: Kynurenine 3-hydroxylase

Application Notes

Optimal dilution of the Kynurenine 3 monooxygenase antibody should be determined by the researcher.

Immunogen

Amino acids 29-42 (KRNFDIDVYEARED) were used as the immunogen for this Kynurenine 3 monooxygenase antibody.

Storage

Aliquot and store the Kynurenine 3 monooxygenase antibody at -20oC.