

Lactadherin Antibody (R36450)

Catalog No.	Formulation	Size
R36450-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q08431
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Lactadherin antibody is available for research use only.



Western blot testing of human MCF7 cell lysate with Lactadherin antibody at 1ug/ml.
Predicted molecular weight: ~43 kDa.

Description

Additional name(s) for this target protein: MFGE8 antibody; Milk Fat Globule-EGF factor 8 protein

Application Notes

Optimal dilution of the Lactadherin antibody should be determined by the researcher.

Immunogen

Amino acids KNNSIPDKQITAS were used as the immunogen for this Lactadherin antibody.

Storage

Aliquot and store the Lactadherin antibody at -20oC.