


LARGE Antibody (R35864)

Catalog No.	Formulation	Size
R35864-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9215
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This LARGE antibody is available for research use only.



Western blot testing of human kidney lysate with LARGE antibody at 1ug/ml. Predicted molecular weight: 88kDa, observed here at ~75kDa.

Description

Additional name(s) for this target protein: Like-acetylglucosaminyltransferase

Application Notes

Optimal dilution of the LARGE antibody should be determined by the researcher.

Immunogen

Amino acids SEADVNSENLQKQ were used as the immunogen for this LARGE antibody.

Storage

Aliquot and store the LARGE antibody at -20oC.