

## LAT2 Antibody (R35762)

Catalog No.	Formulation	Size
R35762-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7462
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This LAT2 antibody is available for research use only.



Western blot testing of human peripheral lymphocytes lysate with LAT2 antibody at 0.01ug/ml. Predicted molecular weight: 27kDa.

## Description

Additional name(s) for this target protein: Linker for activation of T cells family, member 2, Williams-Beuren syndrome chromosomal region 15, WBSCR15

## Application Notes

Optimal dilution of the LAT2 antibody should be determined by the researcher.

## Immunogen

Amino acids RSEKIYQQRSLRED were used as the immunogen for this LAT2 antibody.

## Storage

Aliquot and store the LAT2 antibody at -20oC.