


LCN1 Antibody (R30744)

Catalog No.	Formulation	Size
R30744	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	P31025
Applications	Western Blot : 0.5-1ug/ml IHC (FFPE) : 0.5-1ug/ml
Limitations	This LCN1 antibody is available for research use only.



Western blot testing of LCN1 antibody and Lane 1: Jurkat; 2: COLO320; 3: SCG; 4: HeLa cell lysate

Description

Lipocalin-1, also known as PREALBUMIN, TEAR or VEGP, is a protein that in humans is encoded by the LCN1 gene. Lipocalin-1 is a member of the lipocalin family of small secretory proteins, which is mapped on 9q34.3. It is the primary lipid binding protein in tears and is overproduced in response to multiple stimuli including infection and stress. It may be a marker for chromosome aneuploidy as well as an autoantigen in Sjogren's syndrome. LCN1 inhibits the cysteine-protease papain in vitro, similar to cystatins. And it plays a role in the nonimmunologic defense and in the control of inflammatory processes in oral and ocular tissues. The investigations presented the first clear evidence that LCN1 is induced in infection or inflammation and supported the idea that this lipocalin functions as a physiologic protection factor of epithelia

in vivo.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the LCN1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

An amino acid sequence from the middle region of human LCN1 (HVKDHYIFYCEGELH) was used as the immunogen for this LCN1 antibody.

Storage

After reconstitution, the LCN1 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.