


LIF Antibody (R36251)

Catalog No.	Formulation	Size
R36251-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Predicted Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3976
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This LIF antibody is available for research use only.



Western blot testing of 1) mouse colon and 2) rat colon lysate with LIF antibody at 0.5ug/ml. Observed molecular weight ~23 kDa (unmodified), 32~62 kDa (glycosylated).

Description

Additional name(s) for this target protein: Leukemia inhibitory factor

Application Notes

Optimal dilution of the LIF antibody should be determined by the researcher.

Immunogen

Amino acids TPVNATCAIRHP were used as the immunogen for this LIF antibody.

Storage

Aliquot and store the LIF antibody at -20oC.