


LIPG Antibody (R35611)

Catalog No.	Formulation	Size
R35611-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9388
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This LIPG antibody is available for research use only.



Western blot testing of human liver lysate with LIPG antibody at 0.3ug/ml. Predicted molecular weight: 57kDa but routinely observed at ~70kDa. Both observed bands are blocked by addition of immunizing peptide.

Description

Additional name(s) for this target protein: Endothelial lipase

Application Notes

Optimal dilution of the LIPG antibody should be determined by the researcher.

Immunogen

Amino acids RFNLRTSKDPEHEG were used as the immunogen for this LIPG antibody.

Storage

Aliquot and store the LIPG antibody at -20oC.

References (1)