


LRP5 Antibody (R33982)


Catalog No.	Formulation	Size
R33982-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4041
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 2-6ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This LRP5 antibody is available for research use only.



IHC testing of FFPE human kidney with LRP5 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.



IHC testing of FFPE human liver with LRP5 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.



Western blot testing of human HepG2 cell lysate with LRP5 antibody at 1ug/ml.
Predicted molecular weight ~179 kDa.

Description

Additional name(s) for this target protein: Low density lipoprotein receptor-related protein 5

Application Notes

Optimal dilution of the LRP5 antibody should be determined by the researcher.

Immunogen

Amino acids ERVEKTTGDKRT were used as the immunogen for this LRP5 antibody.

Storage

Aliquot and store the LRP5 antibody at -20oC.