


LSD1 Antibody (R34467)

Catalog No.	Formulation	Size
R34467-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23028
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This LSD1 antibody is available for research use only.



Western blot testing of human kidney lysate with LSD1 antibody at 0.5ug/ml. Expected molecular weight ~110kDa.

Description

Additional name(s) for this target protein: Lysine-specific histone demethylase 1, Amine oxidase (flavin containing) domain 2; AOF2

Application Notes

Optimal dilution of the LSD1 antibody should be determined by the researcher.

Immunogen

Amino acids KHWDQDDDFEFTGSH were used as the immunogen for this LSD1 antibody.

Storage

Aliquot and store the LSD1 antibody at -20oC.