


LTK Antibody / Leukocyte tyrosine kinase receptor [clone 29L03] (RQ8755)

Catalog No.	Formulation	Size
RQ8755	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant RABBIT MONOCLONAL

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	29L03
Purity	Affinity chromatography
UniProt	P29376
Applications	Western Blot : 1:500-1:2000
Limitations	This LTK antibody is available for research use only.



Western blot testing of human Raji cell lysate with LTK antibody. Predicted molecular weight: ~92 kDa, ~86 kDa, ~79 kDa and additional smaller isoforms.

Description

Leukocyte tyrosine kinase receptor is a receptor with a tyrosine-protein kinase activity. Following activation by ALKAL1 or ALKAL2 ligands at the cell surface, transduces an extracellular signal into an intracellular response. [UniProt]

Application Notes

Optimal dilution of the LTK antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human LTK protein was used as the immunogen for the LTK Antibody.

Storage

Store the LTK antibody at -20oC.