


MAG Antibody (R30843)

Catalog No.	Formulation	Size
R30843	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	P20916
Applications	Western Blot : 0.5-1ug/ml (1) Immunohistochemistry (FFPE) : 2-5ug/ml
Limitations	This MAG antibody is available for research use only.



IHC staining of FFPE mouse brain tissue with MAG antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE rat cerebellum tissue with MAG antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.

