


Mboat4 Antibody (R33247)

Catalog No.	Formulation	Size
R33247-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	619373
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Mboat4 antibody is available for research use only.



Western blot testing of mouse pancreas lysate with Mboat4 antibody at 0.5ug/ml. The expected ~50kDa band and the additional ~38kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Membrane bound O-acyltransferase domain containing 4

Application Notes

Optimal dilution of the Mboat4 antibody should be determined by the researcher.

Immunogen

Amino acids CHLGLHYTEYYLGEP were used as the immunogen for this Mboat4 antibody.

Storage

Aliquot and store the Mboat4 antibody at -20oC.