


MBP Antibody (F48188)

Catalog No.	Formulation	Size
F48188-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48188-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P02686
Applications	Western Blot : 1:1000
Limitations	This MBP antibody is available for research use only.



Western blot analysis of MBP antibody and NCI-H460 lysate. Multiple isoforms visualized from 20~37 kDa.

Description

Myelin Basic Protein is a major constituent of the myelin sheath of oligodendrocytes and Schwann cells in the nervous system.

Application Notes

Titration of the MBP antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 177-209 from the human protein was used as the immunogen for this MBP antibody.

Storage

Store at 4oC for up to one month. For long term, aliquot the MBP antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.