


Mcm2 Antibody (F47163)

Catalog No.	Formulation	Size
F47163-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47163-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P97310
Applications	Western Blot : 1:1000
Limitations	This Mcm2 antibody is available for research use only.



Description

Mcm2 acts as a factor that allows the DNA to undergo a single round of replication per cell cycle. Required for the entry in S phase and for cell division.

Application Notes

Titration of the Mcm2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 713-742 from the mouse protein was used as the immunogen for this Mcm2 antibody.

Storage

Aliquot the Mcm2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.