


MDA-5 Antibody / IFIH1 (R35986)

Catalog No.	Formulation	Size
R35986-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64135
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MDA-5 antibody is available for research use only.



IHC testing of FFPE human pancreas with MDA-5 antibody at 10ug/ml. HIER: Microwave with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: Melanoma differentiation associated protein-5, IFIH1, interferon induced with helicase C domain 1

Application Notes

Optimal dilution of the MDA-5 antibody should be determined by the researcher.

Immunogen

Amino acids SNGYSTDENFRYL were used as the immunogen for this MDA-5 antibody.

Storage

Aliquot and store the MDA-5 antibody at -20oC.