


MEFV Antibody / Marenosttrin / Pyrin (R33798)

Catalog No.	Formulation	Size
R33798-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4210
Applications	Immunofluorescence : 10ug/ml
Limitations	This MEFV antibody is available for research use only.



IF/ICC staining of fixed and permeabilized human A431 cells with GADD34 antibody (green) at 10ug/ml, DAPI nuclear stain (blue) and phalloidin (red).

Description

Additional name(s) for this target protein: Mediterranean fever; Marenosttrin; Pyrin

Application Notes

Optimal dilution of the MEFV antibody should be determined by the researcher.

Immunogen

Amino acids EHLKKLRKSGEEQ were used as the immunogen for this MEFV antibody.

Storage

Aliquot and store the MEFV antibody at -20oC.