

Melanoma Associated Antigen Antibody / MAA [clone MAA/1414] (V8899)

Catalog No.	Formulation	Size
V8899-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V8899-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V8899SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2b, kappa
Clone Name	MAA/1414
Purity	Protein A/G affinity
UniProt	Not Known
Localization	Cytoplasmic
Applications	ELISA : 1-5ug/ml (order BSA-free format for coating)
Limitations	This Melanoma Associated Antigen antibody is available for research use only.



Application Notes

Optimal dilution of the Melanoma Associated Antigen antibody should be determined by the researcher.

Immunogen

Proteins extracted from a melanoma were used as the immunogen for the Melanoma Associated Antigen antibody.

Storage

Aliquot the Melanoma Associated Antigen antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.