


MFSD3 Antibody / Major facilitator superfamily domain-containing protein 3 (RQ8703)


Catalog No.	Formulation	Size
RQ8703	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q96ES6
Applications	Western Blot : 1-2ug/ml Flow Cytometry : 1-3ug/million cells ELISA : 0.1-0.5ug/ml
Limitations	This MFSD3 antibody is available for research use only.



Flow cytometry testing of fixed human 293T cells with MFSD3 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= MFSD3 antibody.



Western blot testing of human A549 cell lysate with MFSD3 antibody. Predicted molecular weight ~43 kDa.

Description

MFSD3 belongs to the major facilitator superfamily. The exact functions of MFSD3 remain unknown.

Application Notes

Optimal dilution of the MFSD3 antibody should be determined by the researcher.

Immunogen

An E.coli-derived human recombinant protein (amino acids M1-A381) was used as the immunogen for the MFSD3 antibody.

Storage

After reconstitution, the MFSD3 Antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.