


MIB2 Antibody / Mindbomb 2 (RQ8879)


Catalog No.	Formulation	Size
RQ8879	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q96AX9
Localization	Cytoplasm
Applications	Western Blot : 0.5-1 ug/ml Flow Cytometry : 1-3ug/ml Immunoprecipitation : 2ug per 500ug of lysate
Limitations	This MIB2 antibody is available for research use only.



Flow cytometry testing of fixed and permeabilized human SiHa cells with MIB2 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= MIB2 antibody.



Western blot testing of human 1) A431, 2) HaCaT, 3) SiHa and 4) PC-3 cell lysate with MIB2 antibody. Predicted molecular weight ~110 kDa.



Immunoprecipitation of MIB2 protein from 500ug of human A431 whole cell lysate with 2ug of MIB2 antibody.